

Product Datasheet

GCFC2 Antibody (orb1244273)


Description

GCFC2 Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human C2orf3.
Target	GCFC2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store C2orf3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	C2orf3 antibody can be used for detection of C2orf3 by ELISA at 1:12500. C2orf3 antibody can be used for detection of C2orf3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	89 kDa
Uniprot ID	P16383
NCBI	NP_003194
Dilution Range	C2orf3 antibody can be used for detection of C2orf3 by ELISA at 1:12500. C2orf3 antibody can be used for detection of C2orf3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.




Antibody used in WB on Human Lung at 0.2...

Antibody used in WB on Hum. Fetal Liver ...

Antibody used in WB on Human Jurkat cell...