

Product Datasheet

TCF7L1 Antibody (orb1244294)

Description

TCF7L1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TCF7L1.

Target

TCF7L1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TCF7L1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TCF7L1 antibody can be used for detection of TCF7L1 by ELISA at 1:1562500. TCF7L1 antibody can be used for detection of TCF7L1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

63 kDa

Uniprot ID

[Q9HCS4](#)

NCBI

[NP_112573](#)

Dilution Range


TCF7L1 antibody can be used for detection of TCF7L1 by ELISA at 1:1562500. TCF7L1 antibody can be used for detection of TCF7L1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in WB on Human 293T at 0.2...



Antibody used in WB on Hum. Fetal Liver ...



Antibody used in WB on Hum. Fetal Lung a...