



Product Datasheet

SH3BP2 Antibody (orb1244315)


Description

SH3BP2 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SH3BP2.
Target	SH3BP2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SH3BP2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SH3BP2 antibody can be used for detection of SH3BP2 by ELISA at 1:1562500. SH3BP2 antibody can be used for detection of SH3BP2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	62 kDa
Uniprot ID	P78314
NCBI	NP_003014
Dilution Range	SH3BP2 antibody can be used for detection of SH3BP2 by ELISA at 1:1562500. SH3BP2 antibody can be used for detection of SH3BP2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human kidney at ...



Antibody used in IHC on Human Liver.



Antibody used in WB on Human 721_B cells...