

Product Datasheet

YTHDF3 Antibody (orb1244366)

Description

YTHDF3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human YTHDF3.

Antibody used in WB on Human HeLa at 0.2...

Target

YTHDF3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store YTHDF3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.


Antibody used in WB on Hum. Adult Placen...

Note

For research use only

Application notes

YTHDF3 antibody can be used for detection of YTHDF3 by ELISA at 1:312500. YTHDF3 antibody can be used for detection of YTHDF3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.



Clonality

Polyclonal

MW

64 kDa

Uniprot ID

[Q7Z739](#)

NCBI

[NP_689971](#)

Dilution Range

YTHDF3 antibody can be used for detection of YTHDF3 by ELISA at 1:312500. YTHDF3 antibody can be used for detection of YTHDF3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Antibody used in WB on Hum. Adult Placen...

Expiration Date

12 months from date of receipt.