

## Product Datasheet

### GABARAPL1 Antibody (orb1244395)

## Description

GABARAPL1 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GABARAPL1.

### Target

GABARAPL1

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GABARAPL1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

GABARAPL1 antibody can be used for detection of GABARAPL1 by ELISA at 1:312500. GABARAPL1 antibody can be used for detection of GABARAPL1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

14 kDa

### Uniprot ID

[Q9H0R8](#)

### NCBI


[NP\\_113600](#)

### Dilution Range


GABARAPL1 antibody can be used for detection of GABARAPL1 by ELISA at 1:312500. GABARAPL1 antibody can be used for detection of GABARAPL1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Brain at 5...



Antibody used in IHC on Human Spleen at ...



Antibody used in WB on Human HEK293T cel...