

Product Datasheet

DCAF4 Antibody (orb1244413)

Description

DCAF4 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human WDR21A.

Target

DCAF4

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store WDR21A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

WDR21A antibody can be used for detection of WDR21A by ELISA at 1:312500. WDR21A antibody can be used for detection of WDR21A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

43 kDa

Uniprot ID

[Q8IV10](#)

NCBI

[NP_851937](#)

Dilution Range


WDR21A antibody can be used for detection of WDR21A by ELISA at 1:312500. WDR21A antibody can be used for detection of WDR21A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human Muscle at 0...



Antibody used in WB on Human Jurkat 1 ug...