


Product Datasheet

SNX5 Antibody (orb1244446)

Description


SNX5 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SNX5.
Target	SNX5
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SNX5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SNX5 antibody can be used for detection of SNX5 by ELISA at 1:1562500. SNX5 antibody can be used for detection of SNX5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	47 kDa
Uniprot ID	Q9Y5X3
NCBI	NP_055241
Dilution Range	SNX5 antibody can be used for detection of SNX5 by ELISA at 1:1562500. SNX5 antibody can be used for detection of SNX5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.




See Immunoblot 2 Data and Customer Feedback for more information

Antibody used in WB on Human, Mouse cell...



Antibody used in WB on Human 293T at 0.2...



Antibody used in WB on Human 293T at 1 u...