

Product Datasheet

KPNA4 Antibody (orb1244494)


Description

KPNA4 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KPNA4.
Target	KPNA4
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store KPNA4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	KPNA4 antibody can be used for detection of KPNA4 by ELISA at 1:1562500. KPNA4 antibody can be used for detection of KPNA4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	O00629
NCBI	NP_002259
Dilution Range	KPNA4 antibody can be used for detection of KPNA4 by ELISA at 1:1562500. KPNA4 antibody can be used for detection of KPNA4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HeLa at 0.2...



Antibody used in WB on Mouse Neuro2A at ...



Antibody used in WB on Human MCF7 at 1 u...