

Product Datasheet


GCLM Antibody (orb1244502)

Description

GCLM Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GCLM.
Target	GCLM
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GCLM antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	GCLM antibody can be used for detection of GCLM by ELISA at 1:312500. GCLM antibody can be used for detection of GCLM by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	31 kDa
Uniprot ID	P48507
NCBI	NP_002052
Dilution Range	GCLM antibody can be used for detection of GCLM by ELISA at 1:312500. GCLM antibody can be used for detection of GCLM by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

GCLM




See Immunoblot Data and customer Feedback for more information

Antibody used in WB on Mouse heart at 1:...



Antibody used in WB on Human MCF7 at 0.2...



Antibody used in IHC on Human Liver.