

Product Datasheet

FAM13C Antibody (orb1244547)

Description

FAM13C Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human FAM13C1.

Antibody used in WB on Human Liver at 0....

Target

FAM13C

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store FAM13C1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in WB on Hum. Fetal Muscle...

Note

For research use only

Application notes

FAM13C1 antibody can be used for detection of FAM13C1 by ELISA at 1:312500. FAM13C1 antibody can be used for detection of FAM13C1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Antibody used in WB on Hum. Fetal Brain ...

Clonality

Polyclonal

MW

55 kDa

Uniprot ID

[Q8NE31](#)

NCBI

[NP_001001971](#)

Dilution Range

FAM13C1 antibody can be used for detection of FAM13C1 by ELISA at 1:312500. FAM13C1 antibody can be used for detection of FAM13C1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

16 kDa

40 kDa

31 kDa

22 kDa

FAM13C

90

65

40

31

22

Rabbit Anti-FAM13C1
Sample Type: Human Fetal Muscle
Antibody Concentration: 1 µg/mL

FAM13C

90

65

40

31

22

Rabbit Anti-FAM13C1
Sample Type: Human Fetal Brain
Antibody Concentration: 1 µg/mL