

Product Datasheet


GSR Antibody (orb1244554)

Description

GSR Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GSR.
Target	GSR
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GSR antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	GSR antibody can be used for detection of GSR by ELISA at 1:312500. GSR antibody can be used for detection of GSR by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	P00390
NCBI	NP_000628
Dilution Range	GSR antibody can be used for detection of GSR by ELISA at 1:312500. GSR antibody can be used for detection of GSR by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

Gsr




See Immunoblot 2 Data and customer Feedback for more information

Antibody used in WB on Mouse heart at 1:...



Antibody used in WB on Human HeLa at 0.2...

GSR



Antibody used in WB on zebrafish head at...