

Product Datasheet

ST6GALNAC2 Antibody (orb1244593)

Description

ST6GALNAC2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Heart at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ST6GALNAC2.

Target

ST6GALNAC2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ST6GALNAC2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ST6GALNAC2 antibody can be used for detection of ST6GALNAC2 by ELISA at 1:1562500. ST6GALNAC2 antibody can be used for detection of ST6GALNAC2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

42 kDa

Uniprot ID

[Q9UJ37](#)

NCBI

[NP_006447](#)

Dilution Range

ST6GALNAC2 antibody can be used for detection of ST6GALNAC2 by ELISA at 1:1562500. ST6GALNAC2 antibody can be used for detection of ST6GALNAC2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.