



Product Datasheet

YAP1 antibody (orb12446)


Description	Rabbit polyclonal antibody to YAP1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human YAP1 (NP_001123617.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	75kDa
Uniprot ID	P46937
Entrez	10413
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of PC-12 cell...



Immunofluorescence analysis of NIH-3T3 C...



Immunofluorescence analysis of HeLa cell...