

Product Datasheet

KLB Antibody (orb1244630)

Description

KLB Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KLB.

Target

KLB

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store KLB antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

KLB antibody can be used for detection of KLB by ELISA at 1:12500. KLB antibody can be used for detection of KLB by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

120 kDa

Uniprot ID

[Q86Z14](#)

NCBI


[NP_783864](#)

Dilution Range

KLB antibody can be used for detection of KLB by ELISA at 1:12500. KLB antibody can be used for detection of KLB by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Testis.



Antibody used in WB on Human 721_B at 1 ...



Antibody used in IHC on Human Testis.