

Product Datasheet

PIP5KL1 Antibody (orb1244642)

Description

PIP5KL1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Lung at 0.2...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PIP5KL1.

Target

PIP5KL1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PIP5KL1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PIP5KL1 antibody can be used for detection of PIP5KL1 by ELISA at 1:1562500. PIP5KL1 antibody can be used for detection of PIP5KL1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

22 kDa

Uniprot ID

[Q5T9C9](#)

NCBI

[NP_775763](#)

Dilution Range

PIP5KL1 antibody can be used for detection of PIP5KL1 by ELISA at 1:1562500. PIP5KL1 antibody can be used for detection of PIP5KL1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.