

Product Datasheet

ASB12 Antibody (orb1244751)

Description

ASB12 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Muscle at 0...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ASB12.

Target

ASB12

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ASB12 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ASB12 antibody can be used for detection of ASB12 by ELISA at 1:1562500. ASB12 antibody can be used for detection of ASB12 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

34 kDa

Uniprot ID

[Q8WXX4](#)

NCBI

[NP_569059](#)

Dilution Range

ASB12 antibody can be used for detection of ASB12 by ELISA at 1:1562500. ASB12 antibody can be used for detection of ASB12 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.