

Product Datasheet

CTTN Antibody (orb1244776)

Description

CTTN Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CTTN.

Target

CTTN

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CTTN antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CTTN antibody can be used for detection of CTTN by ELISA at 1:1562500. CTTN antibody can be used for detection of CTTN by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

61 kDa

Uniprot ID

[Q8N707](#)

NCBI

[NP_005222](#)

Dilution Range

CTTN antibody can be used for detection of CTTN by ELISA at 1:1562500. CTTN antibody can be used for detection of CTTN by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.

90kDa
65kDa

40kDa
31kDa
22kDa


Antibody
used in WB
on Human
Brain at 0....

CTTN



Antibody
used in WB
on Hum.
Fetal Lung
a...

CTTN



Antibody
used in WB
on Hum.
Fetal Liver ...