

Product Datasheet

CSNK2A1 Antibody (orb1244779)

Description

CSNK2A1 Antibody

Species/Host

Rabbit

Reactivity

C. elegans, Drosophila, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
721_B cells...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CSNK2A1.

Target

CSNK2A1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CSNK2A1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CSNK2A1 antibody can be used for detection of CSNK2A1 by ELISA at 1:62500. CSNK2A1 antibody can be used for detection of CSNK2A1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

45 kDa

Uniprot ID

[P68400](#)

NCBI

[NP_001886](#)

Dilution Range

CSNK2A1 antibody can be used for detection of CSNK2A1 by ELISA at 1:62500. CSNK2A1 antibody can be used for detection of CSNK2A1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.