

Product Datasheet

ZNF837 Antibody (orb1244828)

Description

ZNF837 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Jurkat 0.2-....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human LOC116412.

Target

ZNF837

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store LOC116412 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

LOC116412 antibody can be used for detection of LOC116412 by ELISA at 1:62500. LOC116412 antibody can be used for detection of LOC116412 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

58 kDa

Uniprot ID

[Q96EG3](#)

NCBI

[XP_375665](#)

Dilution Range

LOC116412 antibody can be used for detection of LOC116412 by ELISA at 1:62500. LOC116412 antibody can be used for detection of LOC116412 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.