

Product Datasheet

LOC388553 Antibody (orb1244829)

Description

LOC388553 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
721_B at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human LOC388553.

Target

LOC388553

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store LOC388553 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

LOC388553 antibody can be used for detection of LOC388553 by ELISA at 1:62500. LOC388553 antibody can be used for detection of LOC388553 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

65 kDa

NCBI

[XP_373809](#)

Dilution Range

LOC388553 antibody can be used for detection of LOC388553 by ELISA at 1:62500. LOC388553 antibody can be used for detection of LOC388553 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.