

Product Datasheet

ZNF705D Antibody (orb1244877)

Description

ZNF705D Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on
Transfected
293T ...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human LOC728957.

Target

ZNF705D

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store LOC728957 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

LOC728957 antibody can be used for detection of LOC728957 by ELISA at 1:62500. LOC728957 antibody can be used for detection of LOC728957 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

35 kDa

Uniprot ID

[A8MUZ8](#)

NCBI

[NP_001034704](#)

Dilution Range

LOC728957 antibody can be used for detection of LOC728957 by ELISA at 1:62500. LOC728957 antibody can be used for detection of LOC728957 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.