



Product Datasheet

MAP1LC3B antibody (orb12449)


Description	Rabbit polyclonal antibody to MAP1LC3B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human LC3B(Q9GZQ8).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa/16kDa
Uniprot ID	Q9GZQ8
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of 293T cell...



Immunohistochemical staining of mouse br...



Immunohistochemical staining of rat brai...