

Product Datasheet

DPCR1 Antibody (orb1244902)

Description

DPCR1 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human DPCR1.
Target	DPCR1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store DPCR1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	DPCR1 antibody can be used for detection of DPCR1 by ELISA at 1:12500. DPCR1 antibody can be used for detection of DPCR1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	55 kDa
Uniprot ID	Q3MIW9
NCBI	CAI18007
Dilution Range	DPCR1 antibody can be used for detection of DPCR1 by ELISA at 1:12500. DPCR1 antibody can be used for detection of DPCR1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

Antibody used in WB on Human HepG2 at 0....