

Product Datasheet

EDAR Antibody (orb1244965)

Description

EDAR Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human EDAR.

Target

EDAR

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store EDAR antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

EDAR antibody can be used for detection of EDAR by ELISA at 1:12500. EDAR antibody can be used for detection of EDAR by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

48 kDa

Uniprot ID

[Q9UNE0](#)

NCBI

[NP_071731](#)

Dilution Range


EDAR antibody can be used for detection of EDAR by ELISA at 1:12500. EDAR antibody can be used for detection of EDAR by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in WB on Human Placenta at...



Antibody used in WB on Hum. Fetal Lung a...



Antibody used in WB on Hum. Fetal Muscle...