

Product Datasheet

SRPRB Antibody (orb1244972)


Description

SRPRB Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SRPRB.
Target	SRPRB
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SRPRB antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SRPRB antibody can be used for detection of SRPRB by ELISA at 1:12500. SRPRB antibody can be used for detection of SRPRB by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	30 kDa
Uniprot ID	Q9Y5M8
NCBI	NP_067026
Dilution Range	SRPRB antibody can be used for detection of SRPRB by ELISA at 1:12500. SRPRB antibody can be used for detection of SRPRB by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HeLa at 0.2...



Antibody used in WB on Human HeLa at 1 u...



Antibody used in WB on Human 721_B at 1 ...