

Product Datasheet

B3GALT1 Antibody (orb1244979)

Description

B3GALT1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human B3GALT1.

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human B3GALT1.

Target

B3GALT1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store B3GALT1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

B3GALT1 antibody can be used for detection of B3GALT1 by ELISA at 1:312500. B3GALT1 antibody can be used for detection of B3GALT1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

38 kDa

Uniprot ID

[Q9Y5Z6](#)

NCBI

[NP_066191](#)

Dilution Range

B3GALT1 antibody can be used for detection of B3GALT1 by ELISA at 1:312500. B3GALT1 antibody can be used for detection of B3GALT1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

80kDa

40kDa

31kDa

22kDa

Antibody used in WB on Human Heart at 0....

B3GALT1

90

65

40

31

22

Rabbit Anti-B3GALT1
Sample Type: Human Fetal Brain
Antibody Concentration: 1µg/mL

Antibody used in WB on Hum. Fetal Brain ...