

Product Datasheet

GALNT16 Antibody (orb1244988)

Description

GALNT16 Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
721_B at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GALNTL1.

Target

GALNT16

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GALNTL1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

GALNTL1 antibody can be used for detection of GALNTL1 by ELISA at 1:62500. GALNTL1 antibody can be used for detection of GALNTL1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

63 kDa

Uniprot ID

[Q8N428](#)

NCBI

[NP_065743](#)

Dilution Range

GALNTL1 antibody can be used for detection of GALNTL1 by ELISA at 1:62500. GALNTL1 antibody can be used for detection of GALNTL1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.