

Product Datasheet

B4GALNT3 Antibody (orb1245027)

Description

B4GALNT3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
HeLa at 0.2...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human B4GALNT3.

Target

B4GALNT3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store B4GALNT3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

B4GALNT3 antibody can be used for detection of B4GALNT3 by ELISA at 1:62500. B4GALNT3 antibody can be used for detection of B4GALNT3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

115 kDa

Uniprot ID

[Q6L9W6](#)

NCBI

[NP_775864](#)

Dilution Range

B4GALNT3 antibody can be used for detection of B4GALNT3 by ELISA at 1:62500. B4GALNT3 antibody can be used for detection of B4GALNT3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.