

Product Datasheet

B3GNT4 Antibody (orb1245041)

Description

B3GNT4 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
293T at 0.2...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human B3GNT4.

Target

B3GNT4

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store B3GNT4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

B3GNT4 antibody can be used for detection of B3GNT4 by ELISA at 1:312500. B3GNT4 antibody can be used for detection of B3GNT4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

42 kDa

Uniprot ID

[Q9C0J1](#)

NCBI

[NP_110392](#)

Dilution Range

B3GNT4 antibody can be used for detection of B3GNT4 by ELISA at 1:312500. B3GNT4 antibody can be used for detection of B3GNT4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.