

## Product Datasheet

### SDF2 Antibody (orb1245064)

## Description

## SDF2 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SDF2.

Antibody used in WB on Human HepG2 at 0....

### Target

SDF2

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SDF2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in WB on Human MCF7 at 1 u...

### Note

For research use only

### Application notes

SDF2 antibody can be used for detection of SDF2 by ELISA at 1:62500. SDF2 antibody can be used for detection of SDF2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

21 kDa

### Uniprot ID

[Q99470](#)

### NCBI

[NP\\_008854](#)

### Dilution Range

SDF2 antibody can be used for detection of SDF2 by ELISA at 1:62500. SDF2 antibody can be used for detection of SDF2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.

