

## Product Datasheet

### POLR3F Antibody (orb1245067)

## Description

POLR3F Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human POLR3F.

Antibody used in WB on Human 721\_B at 0....

### Target

POLR3F

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store POLR3F antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in WB on Human 293T at 1 u...

### Note

For research use only

### Application notes

POLR3F antibody can be used for detection of POLR3F by ELISA at 1:62500. POLR3F antibody can be used for detection of POLR3F by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

36 kDa

### Uniprot ID

[Q9H1D9](#)

### NCBI

[NP\\_006457](#)

### Dilution Range

POLR3F antibody can be used for detection of POLR3F by ELISA at 1:62500. POLR3F antibody can be used for detection of POLR3F by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.

36 kDa

46 kDa

31 kDa

27 kDa

POLR3F

90

65

40

29

22

Protein: Anti-POLR3F  
Sample Type: Human 293T  
Antibody Concentration: 1 µg/mL