


Product Datasheet

GATM Antibody (orb1245080)


Description

GATM Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GATM.
Target	GATM
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GATM antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	GATM antibody can be used for detection of GATM by ELISA at 1:62500. GATM antibody can be used for detection of GATM by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	P50440
NCBI	NP_001473
Dilution Range	GATM antibody can be used for detection of GATM by ELISA at 1:62500. GATM antibody can be used for detection of GATM by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Hum. Fetal Lung a...



Antibody used in WB on Hum. Adult Placen...



Antibody used in WB on Hum. Fetal Heart ...