


Product Datasheet

ACAA1 Antibody (orb1245113)


Description

ACAA1 Antibody


Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ACAA1.
Target	ACAA1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ACAA1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	ACAA1 antibody can be used for detection of ACAA1 by ELISA at 1:1562500. ACAA1 antibody can be used for detection of ACAA1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	P09110
NCBI	NP_001598
Dilution Range	ACAA1 antibody can be used for detection of ACAA1 by ELISA at 1:1562500. ACAA1 antibody can be used for detection of ACAA1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Liver at 1:...



Antibody used in WB on Rat purified pero...



Antibody used in IHC on Human Liver.