

Product Datasheet

SRRD Antibody (orb1245133)

Description

SRRD Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SRRD.

Target

SRRD

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SRRD antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SRRD antibody can be used for detection of SRRD by ELISA at 1:312500. SRRD antibody can be used for detection of SRRD by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

38 kDa

Uniprot ID

[Q9UH36](#)

NCBI

[NP_001013716](#)

Dilution Range

SRRD antibody can be used for detection of SRRD by ELISA at 1:312500. SRRD antibody can be used for detection of SRRD by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

94kDa

41kDa

31kDa

22kDa

Antibody used in WB on Human HepG2 at 0....

SRRD

94

65

40

31

22

Rabbit Anti-SRRD
Sample Type: Human 721_B
Antibody Concentration: 1µg/mL

Antibody used in WB on Human 721_B at 1 ...

SRRD

94

65

40

31

22

Rabbit Anti-SRRD
Sample Type: Human MCF7
Antibody Concentration: 1µg/mL

Antibody used in WB on Human MCF7 at 1 u...