

Product Datasheet

PPP1R13B Antibody (orb1245138)

Description

PPP1R13B Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
HepG2 at
0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PPP1R13B.

Target

PPP1R13B

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PPP1R13B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PPP1R13B antibody can be used for detection of PPP1R13B by ELISA at 1:312500. PPP1R13B antibody can be used for detection of PPP1R13B by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

119 kDa

Uniprot ID

[Q96KQ4](#)

NCBI

[NP_056131](#)

Dilution Range

PPP1R13B antibody can be used for detection of PPP1R13B by ELISA at 1:312500. PPP1R13B antibody can be used for detection of PPP1R13B by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.