

Product Datasheet

PSMC3IP Antibody (orb1245160)

Description

PSMC3IP Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
721_B cells...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PSMC3IP.

Target

PSMC3IP

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PSMC3IP antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PSMC3IP antibody can be used for detection of PSMC3IP by ELISA at 1:1562500. PSMC3IP antibody can be used for detection of PSMC3IP by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

24 kDa

Uniprot ID

[Q9P2W1](#)

NCBI

[NP_037422](#)

Dilution Range

PSMC3IP antibody can be used for detection of PSMC3IP by ELISA at 1:1562500. PSMC3IP antibody can be used for detection of PSMC3IP by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.