

Product Datasheet


CSGALNACT1 Antibody (orb1245175)

Description

CSGALNACT1 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ChGn.
Target	CSGALNACT1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ChGn antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	ChGn antibody can be used for detection of ChGn by ELISA at 1:62500. ChGn antibody can be used for detection of ChGn by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	61 kDa
Uniprot ID	Q8TDX6
NCBI	NP_060841
Dilution Range	ChGn antibody can be used for detection of ChGn by ELISA at 1:62500. ChGn antibody can be used for detection of ChGn by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

Antibody used in WB on Human MCF-7 at 0....



Antibody used in WB on Hum. Fetal Heart ...