

Product Datasheet

TMEM57 Antibody (orb1245182)

Description

TMEM57 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TMEM57.

Target

TMEM57

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TMEM57 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TMEM57 antibody can be used for detection of TMEM57 by ELISA at 1:62500. TMEM57 antibody can be used for detection of TMEM57 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

76 kDa

Uniprot ID

[Q8N5G2](#)

NCBI

[NP_060672](#)


Dilution Range

TMEM57 antibody can be used for detection of TMEM57 by ELISA at 1:62500. TMEM57 antibody can be used for detection of TMEM57 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.


Expiration Date

12 months from date of receipt.

Antibody used in WB on Human MDA-MB-435S...



Antibody used in WB on Hum. Adult Placen...



Antibody used in WB on Hum. Fetal Muscle...