

Product Datasheet

CKAP4 Antibody (orb1245240)


Description

CKAP4 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CKAP4.
Target	CKAP4
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CKAP4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	CKAP4 antibody can be used for detection of CKAP4 by ELISA at 1:312500. CKAP4 antibody can be used for detection of CKAP4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	66 kDa
Uniprot ID	Q07065
NCBI	NP_006816
Dilution Range	CKAP4 antibody can be used for detection of CKAP4 by ELISA at 1:312500. CKAP4 antibody can be used for detection of CKAP4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HeLa at 0.2...



Antibody used in IHC on Human Breast.