

Product Datasheet

ERLIN1 Antibody (orb1245249)

Description

ERLIN1 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ERLIN1.
Target	ERLIN1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ERLIN1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	ERLIN1 antibody can be used for detection of ERLIN1 by ELISA at 1:62500. ERLIN1 antibody can be used for detection of ERLIN1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	O75477
NCBI	NP_006450
Dilution Range	ERLIN1 antibody can be used for detection of ERLIN1 by ELISA at 1:62500. ERLIN1 antibody can be used for detection of ERLIN1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

ERLIN1




Antibody used in WB on Human MCF7 at 1 u...

ERLIN1



Antibody used in WB on Human 721_B at 1 ...



Antibody used in WB on Human HepG2 cells...