

Product Datasheet

GOLGA5 Antibody (orb1245259)

Description

GOLGA5 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GOLGA5.

Antibody used in WB on Human 721_B at 0....

Target

GOLGA5

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GOLGA5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in WB on Human 293T at 1 u...

Note

For research use only

Application notes

GOLGA5 antibody can be used for detection of GOLGA5 by ELISA at 1:312500. GOLGA5 antibody can be used for detection of GOLGA5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

83 kDa

Uniprot ID

[Q8TBA6](#)

NCBI

[NP_005104](#)

Dilution Range

GOLGA5 antibody can be used for detection of GOLGA5 by ELISA at 1:312500. GOLGA5 antibody can be used for detection of GOLGA5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

