

Product Datasheet

LINC00114 Antibody (orb1245331)

Description

LINC00114 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Heart at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NCRNA00114.

Target

LINC00114

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store NCRNA00114 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

NCRNA00114 antibody can be used for detection of NCRNA00114 by ELISA at 1:12500. NCRNA00114 antibody can be used for detection of NCRNA00114 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

15 kDa

Uniprot ID

[Q6XXX2](#)

NCBI

[AAP41721](#)

Dilution Range

NCRNA00114 antibody can be used for detection of NCRNA00114 by ELISA at 1:12500. NCRNA00114 antibody can be used for detection of NCRNA00114 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.