

Product Datasheet

PLD3 Antibody (orb1245398)

Description

PLD3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Heart at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PLD3.

Target

PLD3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PLD3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody
used in WB
on Hum.
Fetal Heart ...

Note

For research use only

Application notes

PLD3 antibody can be used for detection of PLD3 by ELISA at 1:62500. PLD3 antibody can be used for detection of PLD3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

55 kDa

Uniprot ID

[Q8IV08](#)

NCBI

[NP_001026866](#)

Dilution Range

PLD3 antibody can be used for detection of PLD3 by ELISA at 1:62500. PLD3 antibody can be used for detection of PLD3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

91kDa
81kDa

41kDa

31kDa

21kDa

PLD3

90-

85-

40-

35-

22-

Rabbit Anti-PLD3
Sample Type: Human Fetal Heart
Antibody Concentration: 1µg/mL