

Product Datasheet

OR6C75 Antibody (orb1245433)

Description

OR6C75 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Placenta at...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human OR6C75.

Target

OR6C75

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store OR6C75 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

OR6C75 antibody can be used for detection of OR6C75 by ELISA at 1:312500. OR6C75 antibody can be used for detection of OR6C75 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

35 kDa

Uniprot ID

[A6NL08](#)

NCBI

[NP_001005497](#)

Dilution Range

OR6C75 antibody can be used for detection of OR6C75 by ELISA at 1:312500. OR6C75 antibody can be used for detection of OR6C75 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.