

Product Datasheet

EXT2 Antibody (orb1245459)

Description

EXT2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human EXT2.

Target

EXT2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store EXT2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

EXT2 antibody can be used for detection of EXT2 by ELISA at 1:312500. EXT2 antibody can be used for detection of EXT2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

82 kDa

Uniprot ID

[Q93063](#)

NCBI

[NP_000392](#)

Dilution Range

EXT2 antibody can be used for detection of EXT2 by ELISA at 1:312500. EXT2 antibody can be used for detection of EXT2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90 kDa

40 kDa

31 kDa

22 kDa

Antibody used in WB on Human 721_B at 0....

EXT2

90

65

40

28

22

Rabbit Anti-EXT2
Sample Type: Human Hela
Antibody Concentration: 1 µg/mL

Antibody used in WB on Human Hela at 1 u...