

Product Datasheet

SLC33A1 Antibody (orb1245501)

Description

SLC33A1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC33A1.

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC33A1.

Target

SLC33A1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC33A1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLC33A1 antibody can be used for detection of SLC33A1 by ELISA at 1:62500. SLC33A1 antibody can be used for detection of SLC33A1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

61 kDa

Uniprot ID

[O00400](#)

NCBI

[NP_004724](#)

Dilution Range

SLC33A1 antibody can be used for detection of SLC33A1 by ELISA at 1:62500. SLC33A1 antibody can be used for detection of SLC33A1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90 kDa

40 kDa

31 kDa

22 kDa

Antibody used in WB on Human Placenta at...

SLC33A1

90

65

40

31

22

Rabbit Anti-SLC33A1
Sample Type: Human Fetal Brain
Antibody Concentration: 1 µg/mL

Antibody used in WB on Hum. Fetal Brain ...

SLC33A1

90

65

40

31

22

Rabbit Anti-SLC33A1
Sample Type: Human FTL 3
Antibody Concentration: 1 µg/mL

Antibody used in WB on Human 721_B at 1 ...