

Product Datasheet

SLC2A9 Antibody (orb1245514)

Description

SLC2A9 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human SLC2A9.

Target

SLC2A9

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC2A9 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLC2A9 antibody can be used for detection of SLC2A9 by ELISA at 1:1562500. SLC2A9 antibody can be used for detection of SLC2A9 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

56 kDa

Uniprot ID

[Q9NRM0](#)

NCBI

[NP_001001290](#)

Dilution Range

SLC2A9 antibody can be used for detection of SLC2A9 by ELISA at 1:1562500. SLC2A9 antibody can be used for detection of SLC2A9 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90kDa

40kDa

31kDa

22kDa

Antibody used in WB on Human HepG2 at 0....

SLC2A9

90

65

40

31

22

Rabbit Anti-SLC2A9
Sample Type: Human Fetal Lung
Antibody Concentration: 1µg/mL

Antibody used in WB on Hum. Fetal Lung a...