

Product Datasheet

ABCC9 Antibody (orb1245520)

Description

ABCC9 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ABCC9.

Target

ABCC9

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ABCC9 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ABCC9 antibody can be used for detection of ABCC9 by ELISA at 1:1562500. ABCC9 antibody can be used for detection of ABCC9 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

174 kDa

Uniprot ID

[O60706](#)

NCBI

[NP_005682](#)


Dilution Range

ABCC9 antibody can be used for detection of ABCC9 by ELISA at 1:1562500. ABCC9 antibody can be used for detection of ABCC9 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

ABCC9 / Sur2



[See Immunoblot Data and customer Feedback for more information](#)

Antibody used in WB on Mouse brain at 1:...



Antibody used in WB on Human 721_B cells...