

Product Datasheet

Collagen III antibody (PE) (orb124556)

Description

Rabbit polyclonal antibody against Collagen

Species/Host	Rabbit
Reactivity	Canine, Human, Rabbit
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human Collagen alpha 1(III) chain
Target	COL3A1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Application notes	orb124556 is the PE conjugated form of orb10437
Isotype	IgG
Clonality	Polyclonal
MW	117kDa
Uniprot ID	P02461
Expiration Date	12 months from date of receipt.