

Product Datasheet

MPP5 Antibody (orb1245602)

Description

MPP5 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MPP5.

Target

MPP5

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store MPP5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

MPP5 antibody can be used for detection of MPP5 by ELISA at 1:312500. MPP5 antibody can be used for detection of MPP5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

77 kDa

Uniprot ID

[Q8N3R9](#)

NCBI

[NP_071919](#)

Dilution Range


MPP5 antibody can be used for detection of MPP5 by ELISA at 1:312500. MPP5 antibody can be used for detection of MPP5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in WB on Human Jurkat 0.2-...



Antibody used in WB on Hum. Adult Placen...



Antibody used in WB on Hum. Fetal Lung a...